



Submerged 3 G - KIT (Guide)



Features



Component

No	Part Name	Picture	Diameter (Ø)	Length(mm)	Pro' Code
1	Guide Drill		1.8	29	DIGD1808
2	Initial Drill		2.0	33.3	DIID2808
3	Initial Drill		2.0	35.3	DIID2810
4	Initial Drill		2.0	37.3	DIID2812
5	Twist Drill		2.8	33.3	PFD2808
6	Twist Drill		2.8	35.3	PFD2810
7	Twist Drill		2.8	37.3	PFD2812
8	Twist Drill		3.2	33.3	CSFD3208
9	Twist Drill		3.2	35.3	CSFD3210
10	Twist Drill		3.2	37.3	CSFD3212
11	Twist Drill		3.5	33.3	CSFD3508
12	Twist Drill		3.5	35.3	CSFD3510
13	Twist Drill		3.5	37.3	CSFD3512
14	Twist Drill		3.9	33.3	CSFD3908
15	Twist Drill		3.9	35.3	CSFD3910
16	Twist Drill		3.9	37.3	CSFD3912
17	Twist Drill		4.3	33.3	CSFD4308
18	Twist Drill		4.3	35.3	CSFD4310

Component

No	Part Name	Picture	Diameter (Ø)		Length(mm)	Pro' Code
			Type			
19	Twist Drill		4.3	37.3	CSFD4312	
20	Tissue Punch		4.0	25.5	CTP40	
21	Bone Flattening Drill		4.0	31	BFD40	
22	Mountless Driver (2.1 Double Hex)		Machine	35	DMDM14-21	
23	Mountless Driver (2.1 Double Hex)		Hand	37	DMDH22-21	
24	Mountless Driver (2.5 Hex)		Machine	37	DMDM14-25	
25	Mountless Driver (2.5 Hex)		Hand	37	DMDH22-25	
26	Torx Driver		Hand (Short)	21	TDH13	
27	Torx Driver		Hand (Long)	26	TDH18	
28	Torx Driver		Machine (Short)	22	TDM08	
29	Torx Driver		Machine (Long)	32	TDM18	
30	Anchor Drill		1.3	33.4	DAD13	
31	Anchor Screw		1.8	22.2	DAS	
32	Digital Guide Fix		2.1 Non	11.2	DGF21	
33	Digital Guide Fix		2.5 Non	12.3	DGF25	
34	Torque Wrench		Hand	85	TWI/RWI	